

# X36-1 ELECTRIC EXAMINATION TABLE

Registration code: SKD-A





SAIKANG/Medical

# Product data sheet

# X36-1 Features

# AND HIGHLIGHTS

## Technical parameters

External Size	1960×680mm
Height adjustment	470-910mm
Back plate up/down	95°/40°
<ul> <li>Leg plate up/down</li> </ul>	50°/45°
Net weight	68kg
<ul> <li>Safe working load</li> </ul>	160kg
Rated voltage and frequency	AC220V, 50Hz
Input power	500VA

# Technical configuration

🗹 Adjustable arm support	1 pc
☑ Motor	1 pc
Gas spring	2 pcs
☑ 3'' Castor	4 pcs
🗹 Brake Pedal	4 pcs
☑ Hand Controller	1 pc
🗹 Manual Operator	1 pc









#### Framework

Bed frame surface painted with 1.5mm powder coated steel, seamless welded molding by robot, strong corrosion resistance, environmentally-friendly, and bright surface maintain an overall beautiful appearance. With this structure, the whole bed maintains more stable and makes patients feel more safe.



Manual air levers with two air springs to adjust the back and leg plate angle provide users with a quiet and safe environment during diagnosis and avoid interference.



## Hand controller

Equipped with hook-type hand controller, easy to press and adjust optimal height. The hook on the back is convenient to take down for patients and nurses.



## **Round Covered Edge**

The bed platform was covered in 360°. The soft-touch of surface would be safer for users.

#### **Optional Foot-pedal** Controller

Available to operate the entire bed, which can be adjusted to the appropriate height by lightly stepping on it.



### **Mattress**



### **Braking system**

Equipped with Pedal-controlled 3" castors: 4pcs, brake pedals in 2 functions-locked and universal:4 pcs, convenient to move and fix. 360° spinning wheels polyurethane PU-strong, durable, non-slip, strong in pressure resistance, and good in stability. Move freely when sitting on the bed, saving energy and effort.



### **Matching Pillow**

Neoteric design, narrow on the top and wide on the bottom, can be disassembled as required, easy to use. Its well-shaped and elegant design matches with the modern aesthetic, making it harmoniously integrates into the environment.



## Design

The bed surface with two folds of three sections ergonomic design, can realize elevation adjustment according to the patient's height, achieving personalized design to help diagnosis. Bed and chair dualuse, high load-bearing, durable and stable, fearless of heavy pressure impact, to ensure the safety of use.









Design







Manufacturing



Maintenance





### South America Branch

Saikang Medical Technology, Inc W3400 Corrientes, Argentina TEL: +55 (11) 2758-0438 Rua Universal, Hollywood, São Bernardo do Campo, Brazil TEL: +55 (11) 2758-0438

#### East Africa Branch

Saikang Medical Technology Ltd Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

#### South East Asia Branch **US Branch**

Saikang Medical (Philippines) Company Limited Manila, The Philippines TEL: +63 953 044 7864 180.

Healthward Technology, INC 17800 CASTLETON ST STE CITY OF INDUSTRY CA 91748 TEL: +1 (603)416-3063

### Russia Branch

LLC 'SAIKANG' 454004, Russian Federation, Chelyabinsk, Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

### Mid East Branch

Saikang Medical Supplies (Mid East) LLC Musaffah P.O Box 53200, Abu Dhabi, UAE Tel: +971 585710652



Jiangsu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

## **Q** No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China □ +86 18021231901 · S 24/7 ∞ export@saikangmedical.com

() www.saikangmedical.com

# 



