

SKB041-1



SKB041-1

FEATURES AND HIGHLIGHTS





Technical parameters

External size (LxW):	1900x640 mm
Table height adjustment:	550-850 mm
Back plate adjustment:	0-70° (±10°)

Technical configuration

✓ Crank system	1 set
☑ Back plate air spring	1 set
✓ PP side guardrail	1 set
☑ 6" castors with central locking system	4 pcs
☑ Directional wheel	1 pc
☑ Central locking pedal	2 pcs
☑ Directional pedal	2 pcs
✓ IV pole	1 set
✓ IV pole prevision	2 pcs
☑ Oxygen bottle holder	1 pc
✓ Transfer mattress	1 pc
✓ Strips	2 pcs







Metal material

Baosteel from Fortune Global 500 Panasonic Robotic ensure 360° full smooth welding.



Bed frame&chassis

Both with 30*50mm square tube, more strong and stable.



Painting

11 process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12 mm, brightness 60°. Painting can resist 50 kg impact.



Bed platform&guardrai

Bed platform and guardrail are made of high quality PP material.



Plastic

Environmental plastic can be sterilized in 100°C, tensile upto 30 MPa



Back lift

Back lift can be up to 70 degrees, make patient conformtable when needed



Double side castors

Four double side castors of 6 inch diameter, simultaneously lockable by the step pedal at both sides.
TPR tire no worn out after running 30 km, save anti-winding hard shell, united forming without bolts.



Directional center wheel

Can be engaged to facilitate direction control, enabling users to navigate the stretcher easily in hectic and dynamic environments.



Telescopic IV pole

Plastic combined with stainless steel 4 hooks IV pole (Φ16-19), two IV pole previsions seperated on the head&foot part.



Oxygen bottle holder

Two positions for 2L, ϕ 105mm oxygen bottle, one is under head board, another is on the chassis with belt.



Stainless steel handle

Stainless steel handle, more durable.



Strips

With Strips to prevent patients from possibly sliding off the stretcher. Standard thin mattress, 2cm thick, density: 28kg/m3, color: dark blue.

















Design

Manufacturing

Maintenance



USA

 ${\bf Saikang\,Medical\,Technology,Inc.}$ 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497









www.saikangmedical.com















STOCK CODE: 870098

